

Installing the project

1. Download python
 - a. Directly download it via ...
 - b. Use a dedicated environment using conda (or anaconda etc, see their documentation on how to install the packages)
2. Run `/setup.py` to install the required packages
 - a. Alternatively, you can use pip in the command prompt/terminal
`pip install -r ./requirements.txt`

Configuring SQLite databases

3. Download the following from the HITRAN database and order the files following the example of Table 1;
 - a) line-by-line datasets from <https://hitran.org/lbl/>
 - i. need the .par file
 - ii. rename file as “dd_HITRAN_data.par”, where “dd” are two digits of the molecule ID (e.g. 02, 11, ..)
 - b) molparam.txt <https://hitran.org/docs/iso-meta/>
 - c) partition sum lists for all isotopologues of the molecule from <https://hitran.org/docs/iso-meta/>

Table 1: An example for CH₄ on how to install the files downloaded from HITRAN.

/ch4	/raw	/06_HITRAN_dat.par
/ch4	/raw	/molparam.txt
/ch4	/q32.txt	
/ch4	/q33.txt	
/ch4	/q34.txt	
/ch4	/q35.txt	

4. Run the respective python script found under `/ftir_data_processing/generate_database` to obtain the sqlite database similar to shown in Table 2.

Table 2: Output generating sqlite database.

/ch4	/raw	/parsum.dat
/ch4		/ch4_database.sqlite

Test the code

5. Simulate IR absorption spectra using the scripts found in `/ftir_data_processing/simulate_specra/`

If it can simulate the spectra, then I can fit your measured spectra